

PRODUCT SPECIFICATION SHEET

DESCRIPTION: SYSCON HORMONES LEVEL 1

PRODUCT NUMBER: 3CA063

ANALYTE	UNIT	ANALYTE RANGE	ASSAY METHOD
TSH	μIU/mL	1 - 3	NEXIA
T4	μg/dL	2.5 - 4.5	NEXIA
T3	ng/mL	0.6 - 0.8	NEXIA
Free T4	ng/dL	0.4 - 0.8	NEXIA
Free T3	pg/mL	1.8 - 3.2	NEXIA
AFP	ng/mL	15 - 25	NEXIA
CEA	ng/mL	3 - 7	NEXIA
CA 19-9	U/mL	25 - 55	NEXIA
Ferritin	ng/mL	35 - 65	NEXIA
PSA	ng/mL	3.2 - 4.8	NEXIA
PAP	ng/mL	1.4 - 2.6	NEXIA
LH	mIU/mL	3.5 - 6.5	NEXIA
FSH	mIU/mL	2.5 - 5.5	NEXIA
Prolactin	ng/mL	2. - 5.5	NEXIA
HCG	mIU/mL	3 - 7	NEXIA
β-HCG	mIU/mL	2 - 8	NEXIA
Estriol	pg/mL	120 - 180	NEXIA
Progesterone	ng/mL	1.4 - 2.6	NEXIA
Cortisol	μg/dL	3 - 7	NEXIA
IgE	IU/mL	35 - 65	NEXIA
B2MG	mg/L	0.025 - 0.041	NEXIA
IRI	μU/mL	4 - 16	NEXIA
HGH	ng/mL	2.2 - 3.8	NEXIA
CA 125	U/mL	25 - 55	NEXIA
C Peptide	ng/mL	1.5 - 2.5	NEXIA
CKMB	ng/mL	9 - 11	NEXIA
CA 15-3	U/mL	13 - 17	NEXIA
Testosterone	ng/dL	70 - 130	NEXIA

MICROBIOLOGY: ≤ 100 cfu/mL (Colony Forming Units per mL)

STORAGE: 2-8°C

SHELF LIFE: 24 months from date of manufacture

PRODUCT SPECIFICATION SHEET

DESCRIPTION: SYSCON HORMONES LEVEL 2

PRODUCT NUMBER: 3CB063

ANALYTE	UNIT	ANALYTE RANGE	ASSAY METHOD
TSH	μIU/mL	6 - 10	NEXIA
T4	μg/dL	6 - 10	NEXIA
T3	ng/mL	1.4 - 2.0	NEXIA
Free T4	ng/dL	1 - 2	NEXIA
Free T3	pg/mL	3.8 - 6.2	NEXIA
AFP	ng/mL	80 - 120	NEXIA
CEA	ng/mL	13 - 27	NEXIA
CA 19-9	U/mL	70 - 130	NEXIA
Ferritin	ng/mL	200 - 300	NEXIA
PSA	ng/mL	5 - 15	NEXIA
PAP	ng/mL	3.8 - 6.2	NEXIA
LH	mIU/mL	18 - 32	NEXIA
FSH	mIU/mL	8 - 22	NEXIA
Prolactin	ng/mL	14 - 26	NEXIA
HCG	mIU/mL	80 - 120	NEXIA
β-HCG	mIU/mL	60 - 140	NEXIA
Estriol	pg/mL	250 - 350	NEXIA
Progesterone	ng/mL	6 - 12	NEXIA
Cortisol	μg/dL	15 - 25	NEXIA
IgE	IU/mL	220 - 380	NEXIA
B2MG	mg/L	0.09 - 0.15	NEXIA
IRI	μU/mL	25 - 85	NEXIA
HGH	ng/mL	8 - 12	NEXIA
CA 125	U/mL	70 - 130	NEXIA
C Peptide	ng/mL	3 - 7	NEXIA
CKMB	ng/mL	80 - 120	NEXIA
CA 15-3	U/mL	32 - 48	NEXIA
Testosterone	ng/dL	400 - 600	NEXIA

MICROBIOLOGY: ≤ 100 cfu/mL (Colony Forming Units per mL)

STORAGE: 2-8°C

SHELF LIFE: 24 months from date of manufacture

PRODUCT SPECIFICATION SHEET

DESCRIPTION: SYSCON HORMONES LEVEL 3

PRODUCT NUMBER: 3CC063

ANALYTE	UNIT	ANALYTE RANGE	ASSAY METHOD
TSH	μIU/mL	63 - 83	NEXIA
T4	μg/dL	14 - 20	NEXIA
T3	ng/mL	3.2 - 4.8	NEXIA
Free T4	ng/dL	3.4 - 5.0	NEXIA
Free T3	pg/mL	8 - 14	NEXIA
AFP	ng/mL	220 - 340	NEXIA
CEA	ng/mL	40 - 80	NEXIA
CA 19-9	U/mL	210 - 390	NEXIA
Ferritin	ng/mL	600 - 800	NEXIA
PSA	ng/mL	20 - 40	NEXIA
PAP	ng/mL	22 - 34	NEXIA
LH	mIU/mL	80 - 120	NEXIA
FSH	mIU/mL	40 - 80	NEXIA
Prolactin	ng/mL	70 - 130	NEXIA
HCG	mIU/mL	230 - 330	NEXIA
β-HCG	mIU/mL	210 - 350	NEXIA
Estriol	pg/mL	1300 - 1700	NEXIA
Progesterone	ng/mL	14 - 26	NEXIA
Cortisol	μg/dL	30 - 50	NEXIA
IgE	IU/mL	1700 - 2300	NEXIA
B2MG	mg/L	0.25 - 0.31	NEXIA
IRI	μU/mL	100 - 200	NEXIA
HGH	ng/mL	24 - 36	NEXIA
CA 125	U/mL	350 - 650	NEXIA
C Peptide	ng/mL	10 - 20	NEXIA
CKMB	ng/mL	250 - 350	NEXIA
CA 15-3	U/mL	60 - 100	NEXIA
Testosterone	ng/dL	1000 - 1400	NEXIA

MICROBIOLOGY: ≤ 100 cfu/mL (Colony Forming Units per mL)

STORAGE: 2-8°C

SHELF LIFE: 24 months from date of manufacture